

## Product Datasheet

### POLR2G Antibody (orb1269589)

**Description**

POLR2G Antibody

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

FC, WB

**Applications**
**Immunogen**

This POLR2G antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 126-153 amino acids from the C-terminal region of human POLR2G.

**Target**

POLR2G

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**


19 kDa

**Uniprot ID**
**P62487**
**NCBI**
**P62487**
**Dilution Range**


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

**Expiration Date**

12 months from date of receipt.

**Jurkat**  
72


Western blot analysis in Jurkat cell lin...



Flow cytometric analysis of Jurkat cells...