

Product Datasheet

KRTAP13-3 Antibody (orb1269604)

Description

KRTAP13-3 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This KRTAP13-3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 108-135 amino acids from the C-terminal region of human KRTAP13-3.

Target

KRTAP13-3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

19 kDa

Uniprot ID

[Q3SY46](#)

NCBI


[Q3SY46](#)

Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis in Jurkat cell line...



Flow cytometric analysis of Jurkat cells...