

Product Datasheet

CFAP69 Antibody (orb1269605)

Description

CFAP69 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This C7orf63 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 150-178 amino acids from the N-terminal region of human C7orf63.

Target

CFAP69

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


106 kDa

Uniprot ID
A5D8W1
NCBI
A5D8W1
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

MDA-MB435
250


Western blot analysis
in MDA-MB435 cell ...



C7orf63 Antibody
immunohistochemistry
an...