



## Product Datasheet

### SPON1 Antibody (orb1269617)

| Description         | SPON1 Antibody                                                                                                                                                                                                         |                                                                                     |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Species/Host        | Rabbit                                                                                                                                                                                                                 |  |
| Reactivity          | Mouse                                                                                                                                                                                                                  |                                                                                     |
| Conjugation         | Unconjugated                                                                                                                                                                                                           |                                                                                     |
| Tested Applications | IHC-P, WB                                                                                                                                                                                                              |                                                                                     |
| Immunogen           | This VSGP/F-spondin antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 15-43 amino acids from the N-terminal region of human VSGP/F-spondin.                                 | Western blot analysis in mouse kidney ti...                                         |
| Target              | SPON1                                                                                                                                                                                                                  |  |
| Preservatives       | Supplied in PBS with 0.09% (W/V) sodium azide.                                                                                                                                                                         |                                                                                     |
| Form/Appearance     | Liquid                                                                                                                                                                                                                 |                                                                                     |
| Concentration       | batch dependent                                                                                                                                                                                                        |                                                                                     |
| Storage             | Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures. | VSGP/F-spondin Antibody immunohistochemi...                                         |
| Note                | For research use only                                                                                                                                                                                                  |                                                                                     |
| Application notes   | For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50                                                                                                                                             |                                                                                     |
| Isotype             | Rabbit Ig                                                                                                                                                                                                              |                                                                                     |
| Clonality           | Polyclonal                                                                                                                                                                                                             |                                                                                     |
| MW                  | 91 kDa                                                                                                                                                                                                                 |                                                                                     |
| Uniprot ID          | <a href="#">Q9HCB6</a>                                                                                                                                                                                                 |                                                                                     |
| NCBI                | <a href="#">Q9HCB6</a>                                                                                                                                                                                                 |                                                                                     |
| Dilution Range      | For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50                                                                                                                                             |                                                                                     |
| Expiration Date     | 12 months from date of receipt.                                                                                                                                                                                        |                                                                                     |