

Product Datasheet

POMT1 Antibody (orb1269639)

Description

POMT1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This POMT1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 706-735 amino acids from the C-terminal region of human POMT1.

Target

POMT1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:10~50For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

85 kDa

Uniprot ID
Q9Y6A1
NCBI
Q9Y6A1
Dilution Range


For IHC-P starting dilution is: 1:10~50For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



POMT1 Antibody immunohistochemistry anal...



Western blot analysis in ZR-75-1, K562 c...



Flow cytometric analysis of U251 cells (...)