

Product Datasheet

LETM2 Antibody (orb1269653)

Description

LETM2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This LETM2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 101-129 amino acids from the N-terminal region of human LETM2.

Target

LETM2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW




56 kDa

Uniprot ID
[Q2VYF4](#)
NCBI
[Q2VYF4](#)
Dilution Range

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50 For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.


 Western blot analysis
in NCI-H460 cell l...

 LETM2 Antibody
immunohistochemistry
anal...

 Flow cytometric
analysis of NCI-H460
cel...