

Product Datasheet

COX7A2L Antibody (orb1269668)

Description

COX7A2L Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This COX7A2L antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 37-65 amino acids from the Central region of human COX7A2L.

Target

COX7A2L

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

13 kDa

Uniprot ID

[O14548](#)

NCBI


[O14548](#)

Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot
at 1:1000
dilution +
HepG2 ...



Western blot
analysis in
A549 cell line
...



Western blot
analysis in
mouse NIH-
3T3 c...