


Product Datasheet

MBOAT4 Antibody (orb1269685)


Description

MBOAT4 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	This MBOAT4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 258-287 amino acids from the Central region of human MBOAT4.
Target	MBOAT4
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	Q96T53
NCBI	Q96T53
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western Blot at 1:2000 dilution + mouse ...



Western blot analysis in MCF-7 cell line...



Western blot analysis of MBOAT4 Antibody...