

## Product Datasheet

### PLCXD2 Antibody (orb1269709)

## Description

PLCXD2 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

IHC-P, WB

### Applications

### Immunogen

This PLCXD2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 99-127 amino acids from the Central region of human PLCXD2.

### Target

PLCXD2

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

35 kDa

### Uniprot ID

[Q0VAA5](#)

### NCBI


[Q0VAA5](#)

### Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

### Expiration Date

12 months from date of receipt.



Western blot analysis  
in K562 cell line ...



PLCXD2 Antibody  
immunohistochemistry  
ana...