



Product Datasheet

SNX25 Antibody (orb1269726)

Description	SNX25 Antibody	
Species/Host	Rabbit	 Western blot analysis in Jurkat cell lin...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, WB	
Immunogen	This SNX25 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 399-426 amino acids from the Central region of human SNX25.	
Target	SNX25	 Flow cytometric analysis of Jurkat cells...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	98 kDa	
Uniprot ID	Q9H3E2	
NCBI	Q9H3E2	
Dilution Range	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	