

Product Datasheet

NUDT9 Antibody (orb1269729)

Description

NUDT9 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This NUDT9 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 318-346 amino acids from the C-terminal region of human NUDT9.

Target

NUDT9

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

39 kDa

Uniprot ID

[Q9BW91](#)

NCBI


[Q9BW91](#)

Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot at
1:1000 dilution + HeLa
W...



Western blot analysis
in mouse NIH-3T3 c...



NUDT9 Antibody
immunohistochemistry
anal...