

Product Datasheet

SMPD1 Antibody (orb1269734)

Description

SMPD1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This SMPD1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 391-419 amino acids from the C-terminal region of human SMPD1.

Target

SMPD1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:25For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

70 kDa

Uniprot ID

[P17405](#)

NCBI


[P17405](#)

Dilution Range


For FACS starting dilution is: 1:25For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Overlay histogram showing K562 cells sta...



Western Blot at 1:2000 dilution Lane 1: ...



SMPD1 Antibody immunohistochemistry anal...