

## Product Datasheet

### TAF4B Antibody (orb1269769)

## Description

TAF4B Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

FC, WB

### Applications

### Immunogen

This TAF4B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 647-675 amino acids from the C-terminal region of human TAF4B.

### Target

TAF4B

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

91 kDa

### Uniprot ID

[Q92750](#)

### NCBI


[Q92750](#)

### Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

### Expiration Date

12 months from date of receipt.



Western Blot  
at 1:1000  
dilution +  
Ramos ...



Flow  
cytometric  
analysis of  
K562 cells (...)