

Product Datasheet

SPECC1L Antibody (orb1269778)

Description

SPECC1L Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications
Immunogen

This SPECC1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 705-733 amino acids from the Central region of human SPECC1.

Target

SPECC1L

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000
For IHC-P starting dilution is: 1:10~50
For FACS starting dilution is: 1:10~50
For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

125 kDa

Uniprot ID

[Q69YQ0](#)

NCBI

[Q69YQ0](#)

Dilution Range

For WB starting dilution is: 1:1000
For IHC-P starting dilution is: 1:10~50
For FACS starting dilution is: 1:10~50
For IF starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

HL-60

