



Product Datasheet

KLHL35 Antibody (orb1269782)

Description	KLHL35 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This KLHL35 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 301-328 amino acids from the C-terminal region of human KLHL35.	Western blot analysis in MDA-MB231 cell ...  KLHL35 Antibody immunohistochemistry ana...
Target	KLHL35	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	63 kDa	
Uniprot ID	Q6PF15	
NCBI	Q6PF15	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	