

Product Datasheet

MDK Antibody (orb1269808)

Description

MDK Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

This MDK antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 113-143 amino acids from the C-terminal region of human MDK.

Target

MDK

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

16 kDa

Uniprot ID

[P21741](#)

NCBI


[P21741](#)

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis
in mouse kidney ti...



MDK Antibody
immunohistochemistry
analys...



Confocal
immunofluorescent
analysis of M...