


Product Datasheet

B3GAT3 Antibody (orb1269833)

Description	B3GAT3 Antibody	
Species/Host	Rabbit	
Reactivity	Human	Western blot analysis in CEM cell line l...
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This B3GAT3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 293-321 amino acids from the C-terminal region of human B3GAT3.	
Target	B3GAT3	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	37 kDa	
Uniprot ID	O94766	
NCBI	O94766	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	