

Product Datasheet

PLCL2 Antibody (orb1269852)

Description

PLCL2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications
Immunogen

This PLCL2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1099-1126 amino acids from the C-terminal region of human PLCL2.

Target

PLCL2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

126 kDa

Uniprot ID
[Q9UPR0](#)
NCBI
[Q9UPR0](#)
Dilution Range



For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

NCI-H460

250


 Western blot analysis
in NCI-H460 cell l...

 PLCL2 Antibody
immunohistochemistry
anal...

 Confocal
immunofluorescent
analysis of P...