

Product Datasheet

OR4K5 Antibody (orb1269855)

Description

OR4K5 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This OR4K5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 284-313 amino acids from the C-terminal region of human OR4K5.

Target

OR4K5

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

36 kDa

Uniprot ID

[Q8NGD3](#)

NCBI

[Q8NGD3](#)


Dilution Range

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.

HepG2



Western blot analysis in HepG2 cell line...



Flow cytometric analysis of HepG2 cells ...