

Product Datasheet

DDIT3 Antibody (orb1269882)

Description

DDIT3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

FC, IF, WB

Immunogen

This DDIT3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 120-149 amino acids from the C-terminal region of human DDIT3.

Target

DDIT3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:25For WB starting dilution is: 1:2000For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

19 kDa

Uniprot ID

[P35638](#)

NCBI


[P35638](#)

Dilution Range


For FACS starting dilution is: 1:25For WB starting dilution is: 1:2000For IF starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Overlay histogram showing HeLa cells sta...



Overlay histogram showing HeLa cells sta...



Western Blot at 1:2000 dilution + human ...