

Product Datasheet

OR2T3 Antibody (orb1269903)

Description

OR2T3 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

This OR2T3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 280-308 amino acids from the C-terminal region of human OR2T3.

Target

OR2T3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:10~50For WB starting dilution is: 1:1000For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

36 kDa

Uniprot ID

[Q8NH03](#)

NCBI



[Q8NH03](#)

Dilution Range


For IHC-P starting dilution is: 1:10~50For WB starting dilution is: 1:1000For IF starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



OR2T3 Antibody
immunohistochemistry
anal...



Western blot analysis
in HepG2, MDA-
MB45...



Confocal
immunofluorescent
analysis of O...