

Product Datasheet

PABPN1L Antibody (orb1269908)

Description

PABPN1L Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications

Immunogen

This PABPN1L antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 166-194 amino acids from the Central region of human PABPN1L.

Target

PABPN1L

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:500-1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

30 kDa

Uniprot ID

[A6NDY0](#)

NCBI


[A6NDY0](#)

Dilution Range


For WB starting dilution is: 1:500-1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot at 1:500-1:1000 dilution La...



PABPN1L Antibody immunohistochemistry an...



Confocal immunofluorescent analysis of P...