

Product Datasheet

SEMA3G Antibody (orb1269927)

Description

SEMA3G Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This SEMA3G antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 216-244 amino acids from the Central region of human SEMA3G.

Target

SEMA3G

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

87 kDa


Uniprot ID
Q9NS98
NCBI
Q9NS98
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100


Expiration Date

12 months from date of receipt.

MDA-MB231



Western blot analysis
in MDA-MB231 cell ...



SEMA3G Antibody
immunohistochemistry
ana...