




Product Datasheet

TBC1D13 Antibody (orb1269928)

Description	TBC1D13 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	IF, IHC-P, WB	Western blot analysis in CEM cell line l...
Immunogen	This TBC1D13 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 76-105 amino acids from the Central region of human TBC1D13.	
Target	TBC1D13	Western blot analysis in mouse spleen ti...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	TBC1D13 Antibody immunohistochemistry an...
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For IF starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	47 kDa	
Uniprot ID	Q9NVG8	
NCBI	Q9NVG8	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For IF starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	