

Product Datasheet

CBLN2 Antibody (orb1269952)

Description

CBLN2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This CBLN2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 179-207 amino acids from the C-terminal region of human CBLN2.

Target

CBLN2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:10~50For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

24 kDa

Uniprot ID

[Q8IUK8](#)

NCBI


[Q8IUK8](#)

Dilution Range


For IHC-P starting dilution is: 1:10~50For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



CBLN2 Antibody
immunohistochemistry
anal...



Western blot analysis
in mouse brain tis...