

Product Datasheet

PLEKHA4 Antibody (orb1269985)

Description

PLEKHA4 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This PLEKHA4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 21-50 amino acids from the N-terminal region of human PLEKHA4.

Target

PLEKHA4

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

85 kDa

Uniprot ID
[Q9H4M7](#)
NCBI
[Q9H4M7](#)
Dilution Range




For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IHC-P starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

K562

95


 Western blot analysis
in K562 cell line ...

 Flow cytometric
analysis of K562 cells
(...

 PLEKHA4 Antibody
immunohistochemistry
an...