

Product Datasheet

TECRL Antibody (orb1270002)

Description
TECRL Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This SRD5A2L2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 325-354 amino acids from the C-terminal region of human SRD5A2L2.

Target

TECRL

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


42 kDa

Uniprot ID
Q5HYJ1
NCBI
Q5HYJ1
Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

NCI-H460


Western blot analysis in NCI-H460 cell l...



Flow cytometric analysis of NCI-H460 cel...