

## Product Datasheet

### ATG16L2 Antibody (orb1270026)

## Description

ATG16L2 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

IHC-P, WB

### Applications

### Immunogen

This ATG16L2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 277-304 amino acids from the Central region of human ATG16L2.

### Target

ATG16L2

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

69 kDa

### Uniprot ID

[Q8NAA4](#)

### NCBI


[Q8NAA4](#)

### Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

### Expiration Date

12 months from date of receipt.



Western Blot at  
1:1000 dilution Lane  
1: ...



AT  
immunohistochemistry  
analysis in form...