

## Product Datasheet

### HIDE1 Antibody (orb1270074)

## Description

HIDE1 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

WB

### Applications

### Immunogen

This C19orf38 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 129-158 amino acids from the Central region of human C19orf38.

### Target

HIDE1

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

25 kDa

### Uniprot ID

[A8MVS5](#)

### NCBI


[A8MVS5](#)

### Dilution Range

For WB starting dilution is: 1:1000

### Expiration Date

12 months from date of receipt.



Western Blot  
at 1:1000  
dilution +  
293 wh...