

Product Datasheet

PIK3R5 Antibody (orb1270109)

Description
PIK3R5 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This PIK3R5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 763-792 amino acids from the C-terminal region of human PIK3R5.

Target

PIK3R5

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW



97 kDa

Uniprot ID
Q8WYR1
NCBI
Q8WYR1
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis in WiDr cell line ...



PIK3R5 Antibody immunohistochemistry ana...



Flow cytometric analysis of WiDr cells (...)