

Product Datasheet

CCDC54 Antibody (orb1270148)

Description

CCDC54 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This CCDC54 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 91-119 amino acids from the Central region of human CCDC54.

Target

CCDC54

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW



38 kDa

Uniprot ID
Q8NELO
NCBI
Q8NELO
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis
in NCI-H460 cell l...



CCDC54 Antibody
immunohistochemistry
ana...



Flow cytometric
analysis of NCI-H460
cel...