



Product Datasheet

PTCHD4 Antibody (orb1270190)

Description	PTCHD4 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	Western blot analysis in Jurkat cell lin...
Immunogen	This C6orf138 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 563-591 amino acids from the Central region of human C6orf138.	
Target	PTCHD4	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	C6orf138 Antibody immunohistochemistry a...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	96 kDa	
Uniprot ID	Q6ZW05	
NCBI	Q6ZW05	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	