

## Product Datasheet

### LRRC46 Antibody (orb1270211)

**Description**

LRRC46 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

IHC-P, WB

**Applications**
**Immunogen**

This LRRC46 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 38-67 amino acids from the N-terminal region of human LRRC46.

**Target**

LRRC46

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**

35 kDa


**Uniprot ID**
[Q96FV0](#)
**NCBI**
[Q96FV0](#)
**Dilution Range**

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100


**Expiration Date**

12 months from date of receipt.

Y79



Western blot analysis in Y79 cell line I...



LRRC46 Antibody immunohistochemistry ana...