

Product Datasheet

CCDC102A Antibody (orb1270248)

Description

CCDC102A Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Western blot analysis
in mouse bladder t...

Immunogen

This CCDC102A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 401-430 amino acids from the C-terminal region of human CCDC102A.

Target

CCDC102A

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


63 kDa

Uniprot ID
[Q96A19](#)
NCBI
[Q96A19](#)
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.



CCDC102A Antibody
immunohistochemistry
a...