

Product Datasheet

CLEC4F Antibody (orb1270252)

Description

CLEC4F Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This CLEC4F antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 311-339 amino acids from the C-terminal region of human CLEC4F.

Target

CLEC4F

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:500For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

66 kDa

Uniprot ID

[Q8N1N0](#)

NCBI


[Q8N1N0](#)

Dilution Range


For WB starting dilution is: 1:500For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.



Western Blot at 1:500 dilution + CCRF-CE...



CLEC4F Antibody immunohistochemistry ana...