

## Product Datasheet

### ABHD17C Antibody (orb1270256)

**Description**

ABHD17C Antibody

**Species/Host**

Rabbit

**Reactivity**

Mouse

**Conjugation**

Unconjugated

**Tested**

IHC-P, WB

**Applications**
**Immunogen**

This FAM108C1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 92-121 amino acids from the Central region of human FAM108C1.

**Target**

ABHD17C

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000  
For IHC-P starting dilution is: 1:50~100

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**



36 kDa

**Uniprot ID**
[Q6PCB6](#)
**NCBI**
[Q6PCB6](#)
**Dilution Range**


For WB starting dilution is: 1:1000  
For IHC-P starting dilution is: 1:50~100

**Expiration Date**

12 months from date of receipt.



Western blot analysis  
in mouse NIH-3T3 c...



FAM108C1 Antibody  
immunohistochemistry  
a...