

Product Datasheet

TSSK4 Antibody (orb1270280)

Description

TSSK4 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested
Applications

IF, IHC-P, WB

Immunogen

This TSSK4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 226-254 amino acids from the C-terminal region of human TSSK4.

Target

TSSK4

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


37 kDa

Uniprot ID
Q6SA08
NCBI
Q6SA08
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For IF starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of TSSK4 using rab...



TSSK4 Antibody immunohistochemistry anal...



Confocal immunofluorescent analysis of T...