

Product Datasheet

KMT2A Antibody (orb1270287)

Description

KMT2A Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P

Applications

Immunogen

This HRX antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 3879-3908 amino acids from the C-terminal region of human HRX.

Target

KMT2A

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:10~50For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

432 kDa

Uniprot ID

[Q03164](#)

NCBI


[Q03164](#)

Dilution Range


For FACS starting dilution is: 1:10~50For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.



Flow cytometric analysis of CEM cells (r...



Formalin-fixed and paraffin-embedded hum...