

Product Datasheet

GLTSCR1 Antibody (orb1270302)

Description

GLTSCR1 Antibody

MDA-MB435

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This GLTSCR1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1112-1139 amino acids from the Central region of human GLTSCR1.

Target

GLTSCR1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

158 kDa

Uniprot ID

[Q9NZM4](#)

NCBI


[Q9NZM4](#)

Dilution Range

For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis in MDA-MB435 cell ...