

Product Datasheet

CCNG1 Antibody (orb1270305)

Description

CCNG1 Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This CCNG1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 243-272 amino acids from the C-terminal region of human CCNG1.

Target

CCNG1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

34 kDa

Uniprot ID

[P51959](#)

NCBI

[P51959](#)


Dilution Range

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.

m.kidney
95
72



Western blot analysis in mouse kidney ti...



Flow cytometric analysis of HeLa cells (...)