

Product Datasheet

KLF10 Antibody (orb1270308)

Description
KLF10 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This TIEG antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 198-225 amino acids from the Central region of human TIEG.

Target

KLF10

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


53 kDa

Uniprot ID
Q13118
NCBI
Q13118
Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis of TIEG using rabb...



Flow cytometric analysis of 293 cells (r...