



Product Datasheet

THEMIS Antibody (orb1270317)

Description	THEMIS Antibody	
Species/Host	Rabbit	 Western blot analysis in mouse liver tis...
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This THEMIS antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 559-588 amino acids from the C-terminal region of human THEMIS.	 THEMIS Antibody immunohistochemistry ana...
Target	THEMIS	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	73 kDa	
Uniprot ID	Q8N1K5	
NCBI	Q8N1K5	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	