

Product Datasheet

RNF19B Antibody (orb1270321)

Description

RNF19B Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This RNF19B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 369-398 amino acids from the Central region of human RNF19B.

Target

RNF19B

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:500For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

78 kDa

Uniprot ID

[Q6ZMZ0](#)

NCBI


[Q6ZMZ0](#)

Dilution Range


For WB starting dilution is: 1:500For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot
at 1:500
dilution +
THP-1 w...



Flow
cytometric
analysis of
MCF-7 cells
...



Western Blot
at 1:500
dilution +
THP-1 w...