




Product Datasheet

CPSF3L Antibody (orb1270333)

Description	CPSF3L Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IF, IHC-P, WB	Western blot analysis in K562 cell line ...
Immunogen	This CPSF3L antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 91-119 amino acids from the N-terminal region of human CPSF3L.	
Target	CPSF3L	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	CPSF3L Antibody immunohistochemistry ana...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50	Confocal immunofluorescent analysis of C...
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	68 kDa	
Uniprot ID	Q5TA45	
NCBI	Q5TA45	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	