

Product Datasheet

ANKRD39 Antibody (orb1270338)

Description

ANKRD39 Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This ANR39 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 137-165 amino acids from the C-terminal region of human ANR39.

Target

ANKRD39

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

20 kDa

Uniprot ID

[Q53RE8](#)

NCBI

[Q53RE8](#)


Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100


Expiration Date

12 months from date of receipt.

m.cerebellum



Western blot analysis in mouse cerebellu...



ANKRD39 Antibody immunohistochemistry anal...