

Product Datasheet

CUZD1 Antibody (orb1270345)

Description

CUZD1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

This CUZD1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 579-607 amino acids from the C-terminal region of human CUZD1.

Target

CUZD1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

68 kDa

Uniprot ID

[Q86UP6](#)

NCBI


[Q86UP6](#)

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis
in human normal ut...



CUZD1 Antibody
immunohistochemistry
anal...



Confocal
immunofluorescent
analysis of C...