

## Product Datasheet

### HTR2B Antibody (orb1270358)

## Description

## HTR2B Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested

IHC-P, WB

### Applications

### Immunogen

This HTR2B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 408-435 amino acids from the C-terminal region of human HTR2B.

### Target

HTR2B

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

54 kDa

### Uniprot ID

[P41595](#)

### NCBI


[P41595](#)

### Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

### Expiration Date


12 months from date of receipt.



Western Blot at  
1:2000 dilution Lane  
1: ...



Western blot analysis  
in mouse spleen ti...



immunohistochemistry  
analysis in formali...