




## Product Datasheet

### HPR Antibody (orb1270370)

| Description                | HPR Antibody                                                                                                                                                                                                           |                                                                                       |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                 |    |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                                  |                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                           |                                                                                       |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                                                                                                                              |                                                                                       |
| <b>Immunogen</b>           | This HPR antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 212-240 amino acids from the Central region of human HPR.                                                        |                                                                                       |
| <b>Target</b>              | HPR                                                                                                                                                                                                                    |    |
| <b>Preservatives</b>       | Supplied in PBS with 0.09% (W/V) sodium azide.                                                                                                                                                                         |                                                                                       |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                                                                 |                                                                                       |
| <b>Concentration</b>       | batch dependent                                                                                                                                                                                                        |                                                                                       |
| <b>Storage</b>             | Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures. |                                                                                       |
| <b>Note</b>                | For research use only                                                                                                                                                                                                  |  |
| <b>Application notes</b>   | For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:50~100                                                                                                                                            |                                                                                       |
| <b>Isotype</b>             | Rabbit Ig                                                                                                                                                                                                              |                                                                                       |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                             |                                                                                       |
| <b>MW</b>                  | 39 kDa                                                                                                                                                                                                                 |                                                                                       |
| <b>Uniprot ID</b>          | <b>P00739</b>                                                                                                                                                                                                          | Western Blot at 1:8000 dilution + human ...                                           |
| <b>NCBI</b>                | <b>P00739</b>                                                                                                                                                                                                          |                                                                                       |
| <b>Dilution Range</b>      | For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:50~100                                                                                                                                            |                                                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                        |                                                                                       |