


Product Datasheet

IL6ST Antibody (orb1270375)


Description

IL6ST Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	This CD130 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 871-899 amino acids from the C-terminal region of human CD130.
Target	IL6ST
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:25
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	104 kDa
Uniprot ID	P40189
NCBI	P40189
Dilution Range	For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:25
Expiration Date	12 months from date of receipt.



Western Blot
at 1:2000
dilution Lane
1: ...



Western Blot
at 1:2000
dilution Lane
1: ...



Antibody
staining
CD130 in
Human
skeleta...