

Product Datasheet

OR5L2 Antibody (orb1270406)

Description

OR5L2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This OR5L2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 62-90 amino acids from the N-terminal region of human OR5L2.

Target

OR5L2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:50~100For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

35 kDa

Uniprot ID

[Q8NGL0](#)

NCBI


[Q8NGL0](#)

Dilution Range


For IHC-P starting dilution is: 1:50~100For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



OR5L2 antibody immunohistochemistry anal...



Western blot analysis in CEM, HeLa cell ...